





Application Note AN N304

# FT-NIR Spectroscopy for Wet Forage Analysis

The fast analysis of fresh wet forage gains importance for faster screening of the increasing number of samples during harvest. FT-NIR spectroscopy as a rapid and cost effective method, allows determination of multiple nutrients and characteristics of wet and dry forages. The Bruker MPA II with the new large rotating sample cup IN313-S (370 mm diameter) allows easy analysis without sample preparation. The measured sample area is 100 times the area of a static sample and thus enables the analysis of such inhomogeneous samples efficiently and with only one filling.

Alternatively a TANGO or MPA II instrument with a standard 97mm cup can be used with multiple fillings per sample.

## Why scan fresh forage samples?

- Minimal sample preparation
- Faster turnaround time
- More meaningful sample results
- Measure volatile components that aren't driven off in the drying procedure (e.g. ammonia, volatile fatty acids)
- Less chemical changes to sample e.g. Maillard Reaction between protein and sugars

The INGOT calibration packages from Aunir covers various types of forages:

- wet fresh forages and silages
- ensiled grass, silages and maize
- dry roughages (hay, straw) and provides a comprehensive range of parameters:

Dry matter	Oil
D value (organic matter digestibility)	рН
Protein	Lactic acid
Sugar	Volatile fatty acids
Ash	Ammonia
NDF	Starch
ADF	NCGD

# Why choose the Aunir INGOT wet forage calibration package?

- Aunir is the leading independent developer and supplier of NIR calibrations, providing support services to ensure that all calibrations maintain optimum performance.
- INGOT calibrations are built on ISO17025 accredited reference data.
- All data originate from one reference laboratory,
- The INGOT database incorporates many years of seasonal, regional and varietal variation.

#### Hierarchy of INGOT calibration packages for wet forage

Level 4

Wet Forage (Not dried & Not ground)

Silage (Not dried & Not ground)

Fresh Forage (Not dried & Not ground)

Grass Maize Whole Crop Fresh Grass Hay





#### FT-NIR Spectrometers: Bruker Optics offers various FT-NIR spectrometer models for lab, at-line and on-line applications:







MATRIX-F

Process monitoring with probes and measurement heads.

Bruker Optik GmbH

Ettlingen · Deutschland Phone +49 (7243) 504-2000 Fax +49 (7243) 504-2050 info.bopt.de@bruker.com

### **Bruker Optics Inc.**

Billerica, MA · USA Phone +1 (978) 439-9899 Fax +1 (978) 663-9177 info.bopt.us@bruker.com

# Bruker Shanghai Ltd.

Shanghai · China Phone +86 21 51720-890 Fax +86 21 51720-899 info.bopt.cn@bruker.com

#### www.bruker.com/optics

Bruker Optics is continually improving its products and reserves the right to change specifications without notice. © 2018 Bruker Optics BOPT-4001046-02