



TANGO. ANALYSIS TO GO.

• The next generation FT-NIR spectrometer.

Innovation with Integrity

FT-NIR

TANGO CREATES VALUES.Straight forward analyses without delays.

Faster, simpler, more secure - with TANGO your NIR analysis speeds up. TANGO has exactly what users require of an FT-NIR spectrometer suitable for industrial use: robustness, high precision and straightforward operator guidance.

EVERYTHING, WITH EASE.

Fast measurements for high sample throughput, simultaneous evaluation of different components and an intuitive handling, without need for expertise training: TANGO simplifies the NIR analysis and therefore offers maximum security of the results. Specifically targeted for the given application and the demands required of an efficient analyzer. For liquids or for solids, with integrated PC and monitor or as separate analysis station, e.g. for connecting a laptop. Integrated into a network or as a stand alone solution.

ONCE IS ENOUGH. CALIBRATION TRANSFER FOR MORE EFFICIENCY.

For precise analysis, an accurate calibration is crucial. With TANGO, you only need to do it once! All calibrations can be transferred to another spectrometer, also across different system platforms. Easy via LAN or USB stick.



FT-NIR SPECTROMETER

TANGO IS INTUITIVE HANDLING. NIR analysis becomes simpler. And faster. And more secure.

Forget tedious training courses or specialized knowledge. TANGO makes it simple for you: with its intuitive user interface it guides the operators quickly and securely through the work flow. Also untrained staff can operate the spectrometer securely and error free. Good to know: the user interface is available in 17 languages.

into the measuring cup or vial with the measurement being done through glass.

AUTOMATIC: BACKGROUND MEASUREMENTS

Contributing to the high precision of FT-NIR spectrometers is the exact alignment with the background. But TANGO goes one step further: if required, the

A simple touch on the screen starts

ne measurement

states the

The sample is simply filled into the neasurement cup and placed onto the





The product is selected by either pressing

the respective product icon or by using



LOGICAL: SAMPLE HANDLING THE EASY WAY.

When it comes to sample handling, TANGO saves time and prevents operating errors. No sample preparation is necessary, no chemicals, no additional reagents. The sample is simply filled

background measurements are performed automatically without user intervention. This alignment can even take place if a sample is located in the sampling position. The perfect prerequisite for optimal and secure measurements at any time - without (human) errors.

ERGONOMIC: USER COMFORT IN ALL WAYS.

Nowadays, with laboratory space being at a premium, it is difficult to locate a new analyzer in the lab. This isn't an issue with TANGO due to it's small footprint and ability to be positioned in any orientation - it can be operated from all three sides. The monitor swivels in

The progress is displayed by the green status par. After approximately 15 to 30 seconds the easurement is finished

TANGO demonstrates how easy NIR analysis can be, not just during measurements. The housing and touch-screen monitor are rugged and easy to clean - ideal for lab and at-line environments. The exchange of the light source, if required, can be done in a few seconds.

The result view, e.g. as target/actual value mediately

PRACTICAL:







all directions and can even be tilted in various positions to enable an optimum view. Alternatively, it can be placed flat on top of the spectrometer for a perfect viewing from above.

Also the desiccant cartridge, which keeps the inside low on humidity can be changed easily.

ADVANCED EASE OF USE

EFFICIENCY IN EVERY DETAIL.

shows the quality of the analyzed product



One last push of a button: If required, the stored result can be exported via USB at any time.



TANGO IS THE SPECIALIST AT YOUR SIDE.

Innovative technology makes precision a productive factor.



Select the TANGO you need. The reflection analyzer for solids or the transmission analyzer for liquids. With or without integrated PC and monitor.

TRANSMISSION MEASUREMENTS: **HIGHEST PRECISION FOR LIQUIDS.**

TANGO for measuring liquids comes equipped with an active sample heater/ cooler as standard. The temperature of the sample can be selected between +20°C and +80°C. A sensor is permanently checking the temperature of the sample vial. This allows for a fast temperature control of the sample and therefore quick and reliable results.

REFLECTION MEASUREMENTS: ANALYZE SOLIDS EFFICIENTLY.

TANGO for measuring solids contains a gold coated integrating sphere for analysis by diffuse reflection. This sphere ensures highly reproducible measurements of inhomogeneous samples. The diameter of the measurement spot is around 10 mm; the measured sample surface can be enlarged by using sample rotators.

FOR THOSE WHO WANT ADDI-TIONAL RELIABILITY.

TANGO delivers reliable results wherever highest demands are made on raw materials as well as intermediate and finished products and a rapid measurement at low cost must be ensured.

The most important applications are:

- Testing of incoming raw materials
- Process-accompanying analytics
- Final product inspection
- Release analytics

Food and beverage industry

Feed manufacturing

Pharma and Biotechnology





OUTSTANDING FUNCTIONALITY

FOR THOSE WHO PREFER **CUSTOM-MADE APPLICATIONS.**

TANGO offers an efficient and cost-effective analysis for material identification and quantification in a wide variety of application fields. The comprehensive application support of the Bruker specialists helps to select the optimal measuring method as well as appropriate accessories and software solutions.

Typical applications include:

- Food and beverage industry
- Feed manufacturing
- Pharma and biotechnology
- Chemical and petrochemical industry
- Polymer industry

FOR THOSE WHO REQUIRE **COMPLETE SOLUTIONS.**

For markets with standardised products TANGO offers complete solutions with pre-calibrated systems. For the food and feed industry, ready-to-use application packages are available. Please contact your local NIR expert for further information.

Chemical and petrochemical industry Polymer Industry







TANGO IS DIFFERENT.

High-precision measurements independent of operator and environment.



Sturdy housing, precise core. TANGO connects industrial design with modern technology and intuitive handling. The result is measurements with extraordinary precision and reliable stability.

PRECISION SQUARED: CUBE CORNER TECHNOLOGY

Unlike the classical flat mirror technology, Cube Corner mirrors offer a crucial advantage: unbelievable precision of the optical path. Further, the most frequent sources of error such as vibrations do not affect this technology.

LEADING EDGE TECHNOLOGY

When it comes to precision, there is no alternative to Fourier Transform technology. Another advantage is that FT-NIR is independent of ambient light. The sample can be put directly onto the TANGO for measurement without opening any doors or drawers.

OPERATION ON-SITE: PRECISE TECHNOLOGY AND ROBUST DESIGN

Industrial conditions demand a great deal from highly precise measuring instruments. No problem for the TANGO. Due to its special construction it can handle tough conditions and is perfectly suited for rough environments. Besides its encapsulated housing to protect the optics, TANGO is insensitive to vibrations and temperature changes. Behind the design there is a logical philosophy: It is better to exclude errors from the beginning rather than correcting them afterwards.

The easy-to-use software interface is operated via a touch screen panel.

Sample area for solid and liquid samples. The background is recorded automatically if necessary.

2



3

The Internal Validation Unit guarantees that the TANGO is performing at its optimum.

Light source and desiccant cartridge

can be easily exchanged by the

operator.

The Rocksolid[™] interferometer with Cube Corner technology is the heart of the TANGO.

5

Highly reflective gold coated optics ensure optimized light throughput.





TECHNICAL FEATURES

TANGO at a glance:

- Proven FT-NIR technology by Bruker Optics
- Easy-to-use touch screen operation
- Automated background management
- Robust and precise optics
- Small footprint
- User-friendly maintenance
- Direct calibration transfer to and from all Bruker Optics FT-NIR spectrometers
- Data exchange via network
- Complies with 21 CFR Part 11 requirements



All-in-one or stand-alone – TANGO suits your IT structure perfectly.

TANGO IS INNOVATION. Made by Bruker.

Bruker Optics is the leading manufacturer and worldwide supplier of Fourier Transform Infrared, Near Infrared and Raman spectrometers for various industries and applications. For years, we set new standards on the market when it comes to precision and efficiency, ergonomics and ease of operation, consulting and services.



Highest standards in production assure precision and reliability.

HIGHEST QUALITY FROM A **RENOWNED COMPANY: ALWAYS MORE THAN YOU EXPECT**

We are never satisfied with the common market standards. This is where our own research and development departments play a major role: here new ideas are turned into innovative products - in more precision, advanced user comfort and unrivalled reliability.

Plenty of time for personal consultation and customer service guarantee a sustainable and efficient solution.

With service centers all over Europe, North and South America, Asia and Oceania an efficient global technical support is guaranteed. This includes a professional installation as well as qualified and fast after sales service and, if desired remote diagnostics.

MORE THAN STATE-OF-THE-ART TECHNOLOGY: **OUR KNOW-HOW FOR YOU**

Regular training courses and customized local training are vital elements of our know-how transfer. Please contact your local office or look at our web site:

To us, it is obvious that these highest demands are also valid for our production process. High-quality materials, careful workmanship and, if necessary, especially developed production processes and test routines ensure the quality that is common to all Bruker Optics spectrometers. No matter which new products we design, we place the very highest demands on them all. This is why not only the TANGO, but all NIR spectrometers have a 10 years warranty on the moving parts of the interferometer.

WORLDWIDE ON-SITE: WE ARE THERE WHERE YOU NEED US

Bruker's competence is there where our customers need it - from the very first contact. Our application specialists are scientists and engineers which know NIR spectroscopy and spectrometers as well as the customers applications.



KNOW-HOW MEETS SERVICE

www.bruker.com/trainingcourses-optics

A network of international offices and representatives guarantees a competent customer support world wide.



Bruker Optics is ISO 9001 and ISO 13485 certified.

Covered by one or more of the following patents: DE102004025448; DE19940981. Additional patents pending.

www.bruker.com/optics

Bruker Optics Inc.

Billerica, MA · USA Phone +1 (978) 439-9899 Fax +1 (978) 663-9177 info.bopt.us@bruker.com

Bruker Optik GmbH

Ettlingen · Germany Phone +49 (7243) 504-2000 Fax +49 (7243) 504-2050 info.bopt.de@bruker.com Laser class 1 product.

Bruker Shanghai Ltd.

Shanghai · China Phone +86 21 51720-890 Fax +86 21 51720-899 info.bopt.cn@bruker.com

Bruker Optics is continually improving its products and reserves the right to change specifications without notice. © 2017 Bruker Optics BOPT-4000088-06